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Certofix[®] protect

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Literatura

1. Satorius AE, Szafranski J, Pyne D, Ganesan M, Solomon MJ, Newton DW, Bortz DM, Younger JG. Complement c5a generation by staphylococcal biofilms. *Shock*. 2013 Apr; 39(4):336-42
2. Elliott TSJ. The pathogenesis and prevention of intravascular catheter-related infections. In: Hamilton H, Bodenham AR. *Central venous catheters*. Chichester [u.a.]: Wiley-Blackwell 2009; 206-209
3. Ebert T, Smith S, Pancari G, Wu X, Zorman J, Clark D, Cook J, Burns C, Antonello JM, Cope L, Nagy E, Meinke A, McNeely T. Development of a rat central venous catheter model for evaluation of vaccines to prevent *Staphylococcus epidermidis* and *Staphylococcus aureus* early biofilms. *Hum Vaccin*. 2011 Jun; 7(6):630-8
4. Walz JM, Memsoudis SG, Heard SO. Prevention of central venous catheter bloodstream infections. *J Intensive Care Med*. 2010 May-Jun; 25(3):131-8
5. Mermel LA, Farr BM, Sherertz RJ, Raad II, O'Grady N, Harris JS, Craven DE; Infectious Diseases Society of America; American College of Critical Care Medicine; Society for Healthcare Epidemiology of America. Guidelines for the management of intravascular catheter-related infections. *Clin Infect Dis*. 2001 May 1; 32(9):1249-72
6. Yousif A, Jamal MA, Raad I. Biofilm-based central line-associated bloodstream infections. *Adv Exp Med Biol*. 2015; 830:157-79
7. Nakamura I1, Fukushima S2, Hayakawa T2, Sekiya K3, Matsumoto T4. The additional costs of catheter-related bloodstream infections in intensive care units. *Am J Infect Control*. 2015 Oct 1; 43(10):1046-9
8. Frampton GK, Harris P, Cooper K, Cooper T, Cleland J, Jones J, Shepherd J, Clegg A, Graves N, Welch K, Cuthbertson BH. Educational interventions for preventing vascular catheter bloodstream infections in critical care: evidence map, systematic review and economic evaluation. *Health Technol Assess*. 2014 Feb; 18(15):1-365
9. Rosado V, Romanelli RM, Camargos PA. Risk factors and preventive measures for catheter-related bloodstream infections. *J Pediatr (Rio J)*. 2011 Nov-Dec; 87(6):469-77
10. Walder B, Pittet D, Tramèr MR. Prevention of bloodstream infections with central venous catheters treated with anti-infective agents depends on catheter type and insertion time: evidence from a meta-analysis. *Infect Control Hosp Epidemiol*. 2002 Dec; 23(12):748-56
11. The Joint Commission. Preventing Central Line-Associated Bloodstream Infections: A Global Challenge, a Global Perspective. Oak Brook, IL: Joint Commission Resources, May 2012. [Online]
12. Worth LJ, McLaws ML. Is it possible to achieve a target of zero central line associated bloodstream infections? *Curr Opin Infect Dis*. 2012 Dec; 25(6):650-7
13. O'Grady NP, Alexander M, Burns LA et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. CDC. Guidelines for the Prevention of Intravascular Catheter-Related Infections. Centers for Disease Control and Prevention. 2011
14. McLaws ML, Burrell AR. Zero risk for central line-associated bloodstream infection: are we there yet? *Crit Care Med*. 2012 Feb; 40(2):388-93
15. Milstone AM, Reich NG, Advani S, Yuan G, Bryant K, Coffin SE, Huskins WC, Livingston R, Saiman L, Smith PB, Song X. Catheter dwell time and CLABSIs in neonates with PICCs: a multicenter cohort study. *Pediatrics*. 2013 Dec; 132(6):e1609-15
16. Milstone AM, Sengupta A. Do prolonged peripherally inserted central venous catheter dwell times increase the risk of bloodstream infection? *Infect Control Hosp Epidemiol*. 2010 Nov; 31(11):1184-7
17. McLaws ML, Berry G. Nonuniform risk of bloodstream infection with increasing central venous catheter-days. *Infect Control Hosp Epidemiol*. 2005 Aug; 26(8):715-9
18. Guleri A, Kumar A, Morgan RJ, Hartley M, Roberts DH. Anaphylaxis to chlorhexidine-coated central venous catheters: a case series and review of the literature. *Surg Infect (Larchmt)*. 2012 Jun; 13(3):171-4
19. Yasukawa T, Fujita Y, Sari A. Antimicrobial-impregnated central venous catheters. *N Engl J Med*. 1999 Jun 3; 340(22):1762
20. Oda T, Hamasaki J, Kanda N, Mikami K. Anaphylactic shock induced by an antiseptic-coated central venous catheter. *Anesthesiology*. 1997 Nov; 87(5):1242-4
21. Raad I, Hanna H, Jiang Y, Dvorak T, Reitzel R, Chaiban G, Sherertz R, Hachem R. Comparative activities of daptomycin, linezolid, and tigecycline against catheter-related methicillin-resistant *Staphylococcus bacteremic* isolates embedded in biofilm. *Antimicrob Agents Chemother*. 2007 May; 51(5):1656-60
22. Tambe SM, Sampath L, Modak SM. In vitro evaluation of the risk of developing bacterial resistance to antiseptics and antibiotics used in medical devices. *J Antimicrob Chemother* 2001; 47: 589-98
23. Sampath LA, Tambe SM, Modak SM. In vitro and in vivo efficacy of catheters impregnated with antiseptics or antibiotics: evaluation of the risk of bacterial resistance to the antimicrobials in the catheters. *Infect Control Hosp Epidemiol* 2001; 22: 640-6
24. J. Brunke, T. Riemann, I. Roschke, 30 days antimicrobial efficacy of non-leaching central venous catheters (Poster 063), *Critical Care* 2016, Volume 20 Suppl 2

B. Braun Medical s.r.o. | Divize Hospital Care | V Parku 2335/20 | 148 00 Praha 4 | Česká republika
Tel. +420-271 091 111 | Fax +420-271 091 112 | info@bbraun.cz | www.bbraun.cz