
**yuwell**

405A/405B

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
Air Compressive Nebulizer

# User's manual

Please read the user's manual closely!

# Contents

<u>I. Features</u> .....	26
<u>II. Use and maintenance</u> .....	28
<u>III. Others</u> .....	42
<u>IV. EMC Instruction</u> .....	42

# I. Features

## I. Intended use

The 405A/B Air Compressive Nebulizer is a portable medical device for the patient that is intended to aerosolize medications for inhalation. Intended use: Air Compressive Nebulizer is intended for aerosol inhaling therapy in medicinal unit

## II. Product Performance and Index

14. Power:  
Input: AC(100-240)V, (50/60)Hz  
Output: DC12V, 1000mA
15. Maximum nebulization rate:  $\geq 0.2\text{ml/min}$
16. MMD:  $3.9\mu\text{m}\pm 25\%$   
Percentage of nebulized particles with diameter less than  $5\mu\text{m}\geq 65\%$
17. Maximum pressure of compression pump:  $\geq 0.1\text{MPa}$
18. Free flow of compression pump:  $\geq 3.5\text{l/min}$
19. Noise:  $\leq 60\text{dB(A)}$
20. Net weight:  
405A about 260g (unit only)  
405B about 220g (unit only)
21. Size:  
405A Approx. 13cm (l)  $\times$  6.5cm (w)  $\times$  6.8cm (h) 405B Approx. 14.5cm (l)  
 $\times$  8cm (w)  $\times$  5.6cm (h)
22. Contraindication: None
23. Electrical safety requirements: Class II equipment, Type BF applied part
24. Working system: intermittent operation-working for 20 minutes with an interval of 40 minutes
25. Normal working conditions:  
Ambient temperature:  $10^{\circ}\text{C} \sim 40^{\circ}\text{C}$   
Relative humidity: 30%  $\sim$  75%  
Atmospheric pressure: 86kPa  $\sim$  106kPa
26. Transportation and storage environment conditions:  
Ambient temperature:  $-40^{\circ}\text{C} \sim +55^{\circ}\text{C}$   
Relative humidity: 10%  $\sim$  93%, no condensation  
Atmospheric pressure: 70kPa  $\sim$  106kPa

✘ Test conditions: The environmental temperature is 20°C and the humidity is 53%. The test distance from the nebulizing aerosol nozzle is 10cm. (Change with test conditions and test solution)

⚠ Caution:

- ⌚ If the storage temperature is below 5 °C, apply the device only when the temperature meets normal operating conditions
- ⌚ Store this product in a well-ventilated room, and avoid violent vibration in transportation.

### III. Structural Features and Working Principles

4. The pressure source is DC compression pump.
5. Innovative designs, simple operation, compacted size, easy storage and portable.
6. Schematic diagram of working principles

The compressed air generated with the pump inside the main unit, when pumping out from the nozzle, can suck up the medication with the negative pressure between the nozzle and the air tube. The pumped-up medication when impacting on the top of the septum can turn into very fine aerosol before ejecting to the outside.


Figure 1

## IV. Electrical Schematic


Figure 2: Electrical schematic

Electrical maintenance should be operated by professionals.

## II. Use and maintenance

To ensure proper use of this product, please read this instruction manual carefully and operate strictly in accordance with the requirements of the manual. Please contact the supplier or manufacturer if there is any doubt.

### I. Open-package Checking

Users, before installing and commissioning, shall first check whether the appearance of the product is good, and whether the variety and number of attachments is in line with the list of accessories attached to the specifications. If there is any defect, please contact the supplier or manufacturer timely.

### II. Schematic of Each Components and Names of the Whole Machine


Figure 3: Schematic diagram of components (405A)


Figure 4: Schematic diagram of components (405B)


Figure 5: Schematic diagram of the nebulizer kit and accessories


Figure 6: Schematic diagram of the AC adapter

### III. Preparing the Nebulizer for Use

5. Open the cap


6. Add the correct amount of prescribed medication to the medication cup.

⚠ Caution:

- ⚠ Consult with physicians when nebulizing medication.
- ⚠ Pour the correct amount of medication into the medication cup and do not exceed the maximum mark.
- ⚠ The scale on the medication cup is a reference. Please use the scale on your syringe or vial for accurate measurement of medication.
- ⚠ After the medication is added to the medication cup, please properly place the cup and avoid the dumping of the cup and the medication flowing outside.

7. Cover the cap and connect the mouthpiece or mask.


8. Turn the power switch to the off position.

5. Plug the AC adapter with the main unit

6. Plug the adapter into a power outlet


⚠ Caution:

- ⓘ Make sure the AC adapter plug must be fully inserted into the power outlet when in use.
- ⊘ Do not plug or unplug the AC adapter into the electrical outlet with wet hands.

※ Use on holder


Please use the holder when temporarily placing medication cup.

So far, preparation is ready.

#### IV. Steps


9. Turn the power switch to the on position

10. As the main unit starts, connect the air tube to the connector of the main unit and the medication cup respectively.

- 1) Connect one end to the main unit connector.
- 2) Connect the other end to the medication cup connector.


**⚠ Caution:**

① Make sure the air tube is firmly connected to the main unit and the medication cup.

11. Hold the nebulizer kit for treatment.

Inhaling with the mouthpiece or mask

**⚠ Caution:**

① In nebulization, keep the mouthpiece in mouth and close your mouth.


⊘ Do not tilt the nebulizer kit so the angle of the kit is greater than 45 degrees. Medication may flow into the mouth and the spraying may be ineffective.

① If there is anything unusual in using, please stop using it immediately.



ⓘ In nebulization, check whether there is visible aerosol discharging from the medication cup. In the event of irregular spray, please immediately stop the treatment.

ⓘ Be sure to wipe the face after using nebulizer mask, and do not let the medication left in the face.

※ Rotating the cap can properly

12. After using, turn the power switch to the off position.


13. Disconnect the nebulizer kit from the air tube.

14. Check the air tube. No condensation or moisture should remain in the air tube.

15. Disconnect the air tube.

16. Unplug the adapter from the power outlet.

### ⚠Caution:

ⓘ Do not pull the power cord strongly when unplugging it. Please hold the adapter and pull it out.

## V. Cleaning after Each Use

### 1. Cleaning

Following cleaning instructions after each use will prevent any remaining medication in the bottle from drying resulting in the device nebulizing effectively and will help prevent infections.

First, remove the baffle and accessories, discharge remaining medication in medication cup, and do related cleaning or disinfection.

### **1) Baffle, mask and mouthpiece**

Rinse with clean water or immerse in warm water for 15 minutes to ensure hygienic cleaning. Add appropriate amount of vinegar in the water, and then dry naturally in a clean place.

### **2) Main unit and air tube**

Please wipe with moisturized cloth and then dry naturally.

### **3) Methods for removing water inside the air tube:**

- 1) Make sure the air tube is connected to the main unit.
- 2) Remove the air tube from the accessories of medication cup. After the main unit power is turned on, use fingers open and close the air tube ports repeatedly to exhaust water droplets.

#### **⚠Caution:**

- ⓘ Always dispose of any remaining medication in the medication cup after each use. Use fresh medication each time you use the device.
- ⓘ After cleaning, all attachments must be dried and then stored in a clean place.
- ⊘ Do not store the air tube when there is moisture or medication residue inside.

## **2. Disinfection**

Make sure the parts after each use are disinfected. If pollution is serious, please replace parts with new ones timely.

There are two disinfection methods:

### **1) Use a commercially available disinfectant. Please follow the instructions provided by the disinfectant manufacturer.**


- a) Disinfectant soaking should be in accordance with the time required in the instruction manual of disinfectant.

Rinse with clean warm water and allow to air dry in a clean environment.

#### **⚠Caution:**

- ⓘ After the nebulizer kit are removed from the disinfectant, according to instructions on disinfectant, sufficiently wash the units with clean water, and take care to avoid residual disinfectant for disinfectant inhalation may cause worsening of symptoms.
- ⓘ If coloured disinfectants (such as chlorhexidine) are used, after a long time using, there may be different colour changes of units. This is not a problem on the physical properties. When colour change is serious, please purchase separately sold items to replace.

**2) Add enough water in the container, put parts that can be sterilized by boiling, heat up to boiling for 15 to 20 minutes. After boiling, carefully take out the units, shake off excess water and allow to dry in a clean environment.**

**⚠Caution:**

- ⊘ Do not cook the air tube, mask, filter, filter cover and accessories in boiling water or use boiling water to wash to prevent heat deformation.
- ⊘ In the disinfection process, please ensure the medication cup is empty.
- ⊘ When sterilized by boiling, be careful not to burn without water.

Please select the disinfection method listed in the following table:


Part name	Material	Sterilization by boiling	Ingredient: alcohol
			Main disinfectant: disinfection with ethanol
Mouthpiece	PP	○	○
Mask	PVC	✗	○
Nebulizer kit	PP	○	○
Air tube	PVC	✗	✗
Filter	PET	✗	✗

## VI. Replacement of Filter

Typically, a filter sheet shall be replaced every three months.

If the filter has changed colour, replace it with a new one. If medication has spilled on the filter, replace with a new filter immediately.


4. Pull the filter cover to remove from the back side of the compressor

5. Replace the filter.

Remove the old filter sheet with a toothpick or other tools and then install a new one.

6. Install the cover of filter sheet.

### ⚠ Caution:



- ⓘ Before use, make sure the filter sheet is properly installed.
- ⊘ Please use the original filter sheet only and do not use the device if the filter sheet is not installed.
- ⊘ Do not attempt to wash or clean the filter sheet. Damp filters may cause blockages. If water or medication has spilled on the filter, replace with a new filter immediately.
- ⓘ Filter sheet has no positive or negative side.
- ⓘ Please confirm whether the filter sheet is clean and whether there is any stained dust before installing.
- ⓘ To prevent clogs, regularly clean the filter cover sheet with water (non-boiling) or other methods, dry naturally and install properly.

## VII. Maintenance

### ⚠Caution:

- ⊘ Do not use or store when the air tube is bending.
- ⊘ Do not let the main unit, nebulizer kit, AC adapter drop or subjected to strong shocks.
- ⊘ Do not apply force to the air tube and the nebulizer kit to deform.
- ⊘ Do not poke the nebulizer kit with needles or sharp objects.
- ⊘ Continuous use is prevented.
- ⊘ Do not disassemble the main unit or AC adapter or attempt to repair them.
- ⊘ Do not use or store in bathrooms and other damp places.
- ⊘ Do not put the main unit or other parts in extreme heat, cold or direct sunlight.
- ⊘ Do not use benzine, thinner, etc. to clean the air tube.
- ⊕ Please store product away from children, because for it contains small parts that may be swallowed.
- ⊕ Before maintenance, please unplug the AC adapter.

## VIII. Safety Symbols and Meanings in Operating the Device

Symbol	Description	Symbol	Description
	Class II Equipment		Type BF applied part
	Caution		Consult instructions for use
	Recyclable		FRAGILE
	KEEP DRY		KEEP UP
	Manufacturer		Date of manufacturer
	EC-Representative		Serial number
	ON (power, connection to the mains)		OFF (power, disconnection from the mains)
	Medical device		Refer to instructions manual (Background: Blue, Symbol: White)
	This item is compliant with Medical Device Directive 93/42/EEC		

NOTE: IP21 means the limitation degree of protection against access to hazardous parts with a finger and vertically falling water drops.

## IX. Common Failure Analysis and Troubleshooting

No.	Symptom	Analysis	Solution
1	The device is abnormally oud.	The filter cover is incorrectly attached.	Attach the filter cover correctly.
2	The main unit is not operating	Power is not on	Plug the adapter securely into an electrical outlet. Turn the device on.
3	No nebulization when the main unit is operating.	The nebulizer kit has been deformed in disinfection.	Replace the nebulizer kit or air tube.
		The nebulizer kit or air tube has been blocked or deformed	
		Medication runs out.	Add the correct amount of prescribed medication to the medication cup.
4	There are water droplets in air tube.	The medication is excess or the air tube is not dried after cleaning.	After the air tube is connected to the main unit, repeatedly open and close the outlet of the air tube with a finger to get rid of water.
		In nebulization, there is condensation inside the air tube.	

If there are any other questions, please contact the supplier or manufacturer, or please turn to professionals for maintenance.

## X. Notes

### 1. General Safety Precautions:

- 11) Check the main unit and accessories, confirm the existence of any problems and pay particular caution to the following:
  - d) Outlet and the air tube are not damaged.
  - e) Outlet and air tube are not blocked
  - f) The main unit and nebulizer kit can work properly.

- 12) When using this device, there will be some noise and vibration caused by the pump in the compressor. There will also be some noise caused by emission of compressed air from the nebulizer kit. This is normal and does not indicate a malfunction. It will not affect the product life.
- 13) During the use, the main unit may heat up.
- 14) Operate the device only as intended. Do not use the device for any other purposes.
- 15) Do not operate the devices at temperatures greater than +40°C.
- 16) Limit use to 20 minutes at a time, and allow a 40 minutes interval before using the device again.
- 17) Please unplug the AC adapter from the power outlet after use for complete disconnection from the power supply.
- 18) Be careful not to let the power connection socket at the bottom of the main unit and the connection part of AC adapter stain with medication or dust, otherwise, wipe it off quickly with a gauze or cotton swab
- 19) Term of safety use: 3 years (except for wearing consume able parts)
- 20) The caution and warnings in this description are for proper and safe use of the product, and the prevention of harm to the user or others or damages. Warnings and notes are in the following:

Legend	Content
Caution	It indicates that misuse may lead to personal injury or property damage.
	The symbol  indicates mandatory (something that must be followed). Specific mandatory content is in  or is listed in nearby text or drawings. The left picture shows "general mandatory."
	The symbol  indicates prohibition (something that must be followed). Specific prohibited content is in  or is listed in nearby text or drawings. The left picture shows "general prohibition."

-  For patients in coma or sputum patients, please follow the instructions of your physician or licensed healthcare practitioner about the type, dose and regime of medication.
-  Clean and disinfect the nebulizer kit, mask or mouthpiece before using them for the first time after purchase.
-  If the product has not been used for a long period of time, clean and disinfect the nebulizer kit, mask and mouthpiece before using them.

- ⚠ When people use the product, replace the supporting disposable nebulizers (be sure to replace the nebulizer kit and the mouthpiece or mask), or it may cause infection.
- ⊘ When drying the main unit or dried parts, do not use a microwave oven, cutlery dryer or hair dryer.
- ⊘ Do not put fingers or any other objects into the main unit.
- ⊘ Do not use water in the nebulizer for inhaling purposes.
- ⊘ Please use nebulizer kit, air tube and filters manufactured by Yuwell.
- ⊘ Do not use damaged components, mouthpiece, mask and air tube.
- ⊘ Do not cover the main unit and AC adapter with a blanket or towel, or any other type of cover during use.
- ⊘ When the main unit is in normal operation, in addition to cutting off power supply and other necessary operations, do not touch the main unit and AC adapter.
- ⊘ This device is approved for human use only.
- ⊘ Do not use or store the product where the device may be exposed to toxic or volatile gases.
- ⊘ Do not use a mobile phone near the product.
- ⊘ Do not use the product when surrounded by flammable gas.
- ⚠ Provide close supervision when this device is used by children.

## **2. Notes on the power supply**

- ⚠ Wipe off the dust on the AC adapter before use.
- ⊘ Do not use other non-certified AC adapters.
- ⊘ Do not operate the device with a damaged power cord or plug.
- ⊘ Do not wound the cord of AC adapter on the main unit or AC adapter.
- ⚠ When not in use for a long time, unplug the AC adapter.
- ⚠ Unplug the adapter from the main unit after use.
- ⚠ When a power outage occurs during use, turn off the power switch immediately and unplug the AC adapter.
- ⚠ After use and before cleaning, be sure to unplug the AC adapter.
- ⊘ Do not spill water or other liquids on the compressor and AC adapter. These parts are not waterproof. If liquid spills on these parts, immediately unplug the AC adapter and wipe off the liquid with gauze or other soft absorbent material.

## III. Others

### I. List of Accessories

Attachment	Quantity
Main unit	1
Nebulizer kit	1
Ac adapter	1
Filter	5
Instruction manual	1
Warranty	1

If there is any need for parts and accessories, please contact the supplier.

### II. Processing of Waste and Residue

The processing of discarded main unit, accessories and other products should be in accordance with the provisions of the local government.

## IV. EMC Instruction

With the increased number of electronic devices such as PC's and mobile(cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation.

Medical devices should also not interfere with other devices. In order to regulate the requirements for EMC(Electro

Magnetic Compatibility) with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

Medical devices manufactured by our company conform to this EN60601-1-2:2014 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

- ▶ The use of accessories and cables other than those specified by our company, with the exception of cables sold by our company as replacement parts for internal components, may result in increased emission or decreased immunity of the device.
- ▶ The medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device

should be observed to verify normal operation in the configuration in which it will be used.

- ▶ Refer to further guidance below regarding the EMC environment in which the device should be used.

Table1 - Guidance and manufacturer's declaration-electromagnetic emissions.

Guidance and manufacturer's declaration-electromagnetic emissions		
The 405A, 405B Air Compressive Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such environment.		
Emission Test	Compliance	Electromagnetic environment-guidance
Conducted Emission CISPR 11	Group 1 Class B	The 405A, 405B Air Compressive Nebulizer uses RF energy solely for its internal functioning. Therefore, its RF emissions are very low and are not cause interference in proximity of any electronic appliances.  The 405A, 405B Air Compressive Nebulizer is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies building used or domestic purposes.
Radiated Emission CISPR 11	Group 1 Class B	
Harmonic Current Emissions IEC/EN 61000-3-2	Class A	
Voltage Fluctuations and Flicker IEC/EN 61000-3-3	Complies	

Table 2 Guidance and manufacturer's declaration – electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity			
The 405A, 405B Air Compressive Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such environment.			
Immunity Test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment-guidance
Electrostatic Discharge IEC/EN 61000-4-2	±8KV contact ±15KV air	No degradation of function	Floor should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrostatic Fast Transients /Burst	±2 KV 100 KHz	No degradation of function	Mains power quality should be that of a typical commercial environment or hospital environment.

IEC/EN 61000-4-4	Repetition frequency		
Surges IEC/EN 61000-4-5	±1 KV Line to line	No degradation of function	Mains power quality should be that of a typical commercial environment or hospital environment.
Voltage Dips IEC/EN 61000-4-11	0%U <sub>T</sub> , 0.5 cycle at 0°, 45°90°135°180° 225°270° and 315°;  0%U <sub>T</sub> , 1 cycle at 0°; 70% U <sub>T</sub> ,  25(50Hz)/30 (50Hz)cycle at 0°; 0% U <sub>T</sub> ,  250(50Hz)/300 (60Hz)cycle at 0°;	No degradation of function	Mains power quality should be that of a typical commercial environment or hospital environment. If the user of the portable phlegm suction unit 405A, 405B requires continued operation during power mains interruption, it is recommended that the product be powered from an uninterruptible power supply or battery.
Voltage Interruptions IEC/EN 61000-4-11	0%U <sub>T</sub> ; 250/300cycle		Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Power frequency Magnetic Field IEC/EN 61000-4-8	30 A/m	No degradation of function	
Note: U <sub>T</sub> is the A.C. mains voltage prior to application of the test level			

Table 3 Guidance and manufacturer's declaration – electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity			
The 405A, 405B Air Compressive Nebulizer is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such environment.			
Immunity Test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment-guidance
Conducted RF Disturbances IEC/EN 61000-4-6	3V 150 kHz - 80 MHz	3V rms	The portable and mobile RF communication devices, including

<p>Radiated RF Disturbances IEC/EN 61000-4-3</p>	<p>10V/m80MHz- 2.7GHz 80%AMat1kHz</p>	<p>10 V/m</p>	<p>cables, must not be used closer to the 405A, 405B, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance  <math>d=1.2 \sqrt{P}</math>  <math>d=1.2 \sqrt{P}</math> 80 MHz to 800 MHz  <math>d=2.3 \sqrt{P}</math> 800 MHz to 2.7 GHz</p> <p>where P is the maximum nominal output voltage of the transmitter in watts (W) depending on the manufacturer of the transmitter and the recommended separation distance in meters (m). The intensity of the field from the fixed RF transmitters, as determined by an Electromagnetic study of the site), could be lower than the level of conformity of each frequency interval). It is possible to check for interference in proximity to devices identified by the following symbol:</p> 
--	---	---------------	--

**Note1:** At 80 MHz and 800 MHz, the higher frequency is applied.

**Note2:** These guidelines may not be applicable in all situations.

Electromagnetic propagation is influenced by the absorption and by reflection from buildings, objects and people.

a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen. To establish an electromagnetic environment generated by fixed RF transmitters, an electromagnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored.

If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary.


b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.

Table 4 Recommended separation distance

Recommended separation distance between portable and mobile RF communications equipment and the 405A, 405B.			
The 405A, 405B Air Compressive Nebulizer is intended to operate in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of 405A, 405B can help prevent electromagnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and 405A, 405B as recommended below, according to the maximum output power of the communications equipment.			
Maximum nominal output power of the Transmitter(W)	Separation distance according to frequency of transmitter in meter		
	150 kHz~80 MHz $d=1.2 \sqrt{P}$	80 MHz~800 MHz $d=1.2 \sqrt{P}$	800 MHz~2.7 GHz $d=2.3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note1: At 80 MHz and 800 MHz, the higher frequency is applied.			
Note2: These guidelines may not be applicable in all situations.			
Electromagnetic propagation is influenced by the absorption and by reflection from buildings, objects and people.			

The company reserves the rights to change the product technology and appearance. Please don't mind!

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